


# Surgivisio Instrumentation



## ● Features:

- Single-use
- Pre-calibrated in production
- Delivered sterile

### Shouldered broaches for spinous process fixation



Broaches external diameter: Ø 2mm  
Length: 100mm  
Shouldering: Ø 1mm by 4mm

### Patient reference + 2 wedges




### 3D imaging Phantom



### 3D imaging Offset Phantom



### Patient tracker for navigation



### Navigated targeting needles bevel tip / trocar tip




Despite broaches, each instrument is identified by an identification card which can be inserted into the system. The instrument is automatically recognized.



# ● For each surgery:

## Imaging



3D Imaging Set SPX1  
A00018


OR



3D Offset Imaging Set SPX1  
A00038



## Broach Set (x4)



Broach set SPX1(x4)  
A00021



## Navigation



Navigation Set SPX1,  
Targeting Needle Bevel Tip  
- 11G, 125mm  
A00019

OR




Navigation Set SPX1,  
Targeting Needle Trocar Tip  
- 11G, 125mm  
A00029

## Additional Navigated targeting needles



Navigated Targeting Needle  
SPX1, Bevel Tip

+



identification card


**Navigated Targeting Needle  
SPX1, Bevel Tip - 11G, 125mm  
A00020**

**OR**



Navigated Targeting Needle  
SPX1, Trocar Tip

+



identification card

**Navigated Targeting Needle  
SPX1, Trocar Tip - 11G, 125mm  
A00028**

## Protective covers

Remote panel sterile protective cover

C-arm sterile protective cover

Camera handle sterile protective cover

## References:

<b>A00018</b>	3D Imaging Set SPX1(class Is)
<b>A00038</b>	3D Offset Imaging Set SPX1(class Is)
<b>A00021</b>	Broach set SPX1(x4)(class IIa)
<b>A00019</b>	Navigation Set SPX1, Targeting Needle Bevel Tip - 11G, 125mm (class IIa)
<b>A00029</b>	Navigation Set SPX1, Targeting Needle Trocar Tip - 11G, 125mm (class IIa)
<b>A00020</b>	Navigated Targeting Needle SPX1, Bevel Tip - 11G, 125mm (class IIa)
<b>A00028</b>	Navigated Targeting Needle SPX1, Trocar Tip - 11G, 125mm (class IIa)
<b>I00011</b>	Remote panel sterile protective cover (class Is)
<b>I00009</b>	C-arm sterile protective cover (class Is)
<b>I00010</b>	Camera handle sterile protective cover (class Is)

*For more information on this device, please refer to the User Manual (Ref. AA-0711).*

Our Customer Service is available for any order or additional information at  
[service@ecential-robotics.com](mailto:service@ecential-robotics.com)



ECENTIAL ROBOTICS

**ecential**robotics